

Ordering Information

OTW (Over The Wire)

► Fully covered

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TR1603FP	16	3	18	60
TR1604FP		4		
TR1605FP		5		
TR1606FP		6		
TR1607FP		7		
TR1608FP		8		
TR1609FP		9		
TR1610FP		10		
TR1803FP	18	3		
TR1804FP		4		
TR1805FP		5		
TR1806FP		6		
TR1807FP		7		
TR1808FP		8		
TR1809FP		9		
TR1810FP		10		
TR2003FP	20	3		
TR2004FP		4		
TR2005FP		5		
TR2006FP		6		
TR2007FP		7		
TR2008FP		8		
TR2009FP		9		
TR2010FP		10		
TR2203FP	22	3		
TR2204FP		4		
TR2205FP		5		
TR2206FP		6		
TR2207FP		7		
TR2208FP		8		
TR2209FP		9		
TR2210FP		10		
TR2403FP	24	3		
TR2404FP		4		
TR2405FP		5		
TR2406FP		6		
TR2407FP		7		
TR2408FP		8		
TR2409FP		9		
TR2410FP		10		

► Uncovered

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TR1603P	16	3	16	60
TR1604P		4		
TR1605P		5		
TR1606P		6		
TR1607P		7		
TR1608P		8		
TR1609P		9		
TR1610P		10		
TR1803P	18	3		
TR1804P		4		
TR1805P		5		
TR1806P		6		
TR1807P		7		
TR1808P		8		
TR1809P		9		
TR1810P		10		
TR2003P	20	3		
TR2004P		4		
TR2005P		5		
TR2006P		6		
TR2007P		7		
TR2008P		8		
TR2009P		9		
TR2010P		10		

* Coding

TR--P : P-distal side flare

* Coding

TR--FP : F-fully covered / P-distal side flare